

Prepare the MEG System

Prepare MSR and choose gantry position

Start *megacq*

Set Project and add Participant anonymous information

Check MEG sensor quality

Prepare stimulus equipment and software, triggers recorded etc.

Take an empty room recording

Prepare the participant

Explain preparation procedure

Explain experimental protocol

Documentation read and signed

Metal detection for sources of artifact

Apply ECG, EOG electrodes, EEG/EMG as necessary

Check electrode impedances (SIGGI-II multimeter)

Apply HPI coils

Digitise landmarks, HPI coils, headshape, EEG (if used)

Prepare the session

Position participant in gantry chair

Plug in electrodes, EEG, audio etc.

Position table, set up EyeLink (if used)

Confirm instructions with participant

Start *megacq*, measure/inspect HPI

Adjust gantry, EyeLink hardware (as necessary)

“Accept” HPI, check MEG sensor quality

Confirm all OK with participant

Acquire data

Select *chPI*, Record raw to acquire data

Start experiment, stimulation (as necessary)

Monitor participant, experimental responses, and equipment

Monitor MEG sensor quality

Stop the recording, save data to disk

Save any Behavioural data (as necessary)

Finalise the session

Unplug electrodes, EEG, move table, EyeLink hardware

Help participant out of MSR, debrief/answer any queries

Remove electrodes, EEG, Offer hair washing

Copy off data to RDS (if time permits)

Clean Up, Tidy Up